



# Automatic Cleaning Sequence for: Glass Pipette-Washer Type: MPS

The MPS system is applicable to all validation requirements



**Program 1:** cleans the pipettes in a 95°C alkaline bath, followed by two cold rinses plus one hot rinse at 90°C  
**Program 2:** cleans the pipettes in a 70°C alkaline bath, followed by one cold rinse plus one hot rinse 90°C  
**Program 1 and Program 2** may be adjusted to almost any other sequence even at a later time

**H. Hölzel GmbH**  
 Bahnhofstrasse 23, 85457 Hörlkofen near Erding  
 Tel : (0049) 8122 – 9748 - 0, Fax: (0049) 8122 – 9748 - 21,  
 E-mail: [info@hoelzel-gewo.com](mailto:info@hoelzel-gewo.com), Homepage: [www.hoelzel-gewo.com](http://www.hoelzel-gewo.com)